PASS Sample Size Software NCSS.com

Chapter 233

Non-Zero Null Tests for the Difference of Two Proportions in a Cluster-Randomized Design

Introduction

This module provides power analysis and sample size calculation for inequality tests of the difference with a non-zero null hypothesis in two-sample, cluster-randomized designs in which the outcome is binary.

Technical Details

Our formulation comes from Donner and Klar (2000). Denote a binary observation by Y_{gkm} where g=1 or 2 is the group, $k=1,2,...,K_g$ is a cluster within group g, and $m=1,2,...,M_g$ is an individual in cluster k of group g. The results that follow assume an equal number of individuals per cluster. When the number of subjects from cluster to cluster are about the same, the power and sample size values should be fairly accurate. In these cases, the average number of subjects per cluster can be used.

The statistical hypothesis that is tested concerns the difference between the two group proportions, p_1 and p_2 . When necessary, we assume that group 1 is the treatment group and group 2 is the control group. With a simple modification, all of the large-sample sample size formulas that are listed in the module for testing two proportions can be used here.

When the individual subjects are randomly assigned to one of the two groups, the variance of the sample proportion is

$$\sigma_{S,g}^2 = \frac{p_g (1 - p_g)}{n_o}$$

When the randomization is by clusters of subjects, the variance of the sample proportion is

$$\sigma_{C,g}^2 = \frac{p_g (1 - p_g) (1 + (m_g - 1)\rho)}{k_g m_g}$$
$$= \sigma_{S,g}^2 [1 + (m_g - 1)\rho]$$
$$= F_{g,\rho} \sigma_{S,g}^2$$

The factor $\left[1+\left(m_g-1\right)\rho\right]$ is called the *inflation factor*. The Greek letter ρ is used to represent the *intracluster correlation coefficient (ICC)*. This correlation may be thought of as the simple correlation between any two subjects within the same cluster. If we stipulate that ρ is positive, it may also be interpreted as the proportion of total variability that is attributable to differences between clusters. This value is critical to the sample size calculation.

The asymptotic formulas that were used in comparing two proportions (see Chapter 205, "Tests for the Difference Between Two Proportions with a Non-Zero Null Hypothesis") may be used with cluster-randomized designs as well, as long as an adjustment is made for the inflation factor.

Power Calculations

A large sample approximation may be used that is most accurate when the values of n_1 and n_2 are large. The large approximation is made by replacing the values of \hat{p}_1 and \hat{p}_2 in the z statistic with the corresponding values of p_1 and p_2 under the alternative hypothesis, and then computing the results based on the normal distribution.

Note that in this case, exact calculations are not possible.

Procedure Options

This section describes the options that are specific to this procedure. These are located on the Design and Options tabs. For more information about the options of other tabs, go to the Procedure Window chapter.

Design Tab

The Design tab contains the parameters associated with this test such as the proportions, sample sizes, alpha, and power.

Solve For

Solve For

This option specifies the parameter to be solved for using the other parameters. The parameters that may be selected are *Power*, *Sample Size* (*K1*), *Sample Size* (*M1*), *Effect Size*, and *ICC*. Under most situations, you will select either *Power* or *Sample Size* (*K1*).

Select *Sample Size* (*K1*) when you want to calculate the number of clusters per group needed to achieve a given power and alpha level.

Select *Power* when you want to calculate the power of an experiment.

Test

Alternative Hypothesis

Specify the direction of the alternative hypothesis.

- One-Sided (H1: P1 P2 > D0)
 Refers to an upper-tailed, one-sided test in which the alternative hypothesis is of the type H1: P1 P2 > D0.
- One-Sided (H1: P1 P2 < D0)

 Refers to a lower-tailed, one-sided test in which the alternative hypothesis is of the type H1: P1 P2 < D0.

• Two-Sided (H1: P1 - P2 ≠ D0)

Refers to a two-sided test in which the alternative hypothesis is of the type H1: P1 - P2 \neq D0.

For one-sided tests, the direction you select must match the values entered for P1.0 (or D0), P1.1 (or D1), and P2 you enter. For example, if you select H1: P1 - P2 < D0, then the value(s) for D1 must be less than the value(s) for D0 (and the value(s) for P1.1 must be less than P1.0).

Test Type

Specify which test statistic is used in searching and reporting. We recommend the likelihood score test.

Power and Alpha

Power

This option specifies one or more values for power. Power is the probability of rejecting a false null hypothesis, and is equal to one minus Beta. Beta is the probability of a type-II error, which occurs when a false null hypothesis is not rejected.

Values must be between zero and one. Historically, the value of 0.80 (Beta = 0.20) was used for power. Now, 0.90 (Beta = 0.10) is also commonly used.

A single value may be entered here or a range of values such as 0.8 to 0.95 by 0.05 may be entered.

Alpha

This option specifies one or more values for the probability of a type-I error. A type-I error occurs when a true null hypothesis is rejected.

Values must be between zero and one. Historically, the value of 0.05 has been used for alpha. This means that about one test in twenty will falsely reject the null hypothesis. You should pick a value for alpha that represents the risk of a type-I error you are willing to take in your experimental situation.

You may enter a range of values such as 0.01 0.05 0.10 or 0.01 to 0.10 by 0.01.

Sample Size – Group 1 (Treatment)

K1 (Clusters in Group 1)

Enter a value (or range of values) for the number of clusters in group one. You may enter a range of values such as 10 to 20 by 2. The sample size for this group is equal to the number of clusters times the number of subjects per cluster.

M1 (Items per Cluster in Group 1)

This is the average number of items (subjects) per cluster in group one. This value must be a positive number that is at least 1. You can use a list of values such as 100 150 200.

Sample Size - Group 2 (Reference)

K2 (Clusters in Group 2)

This is the number of clusters in group two. The sample size for this group is equal to the number of clusters times the number of subjects per cluster. This value must be a positive number.

If you simply want a multiple of the value for group one, you would enter the multiple followed by K1, with no blanks. If you want to use K1 directly, you do not have to pre-multiply by 1. For example, all of the following are valid entries: $10 \ K1 \ 2K1 \ 0.5K1$.

You can use a list of values such as 10 20 30 or K1 2K1 3K1.

M2 (Items per Cluster in Group 2)

This is the number of items (subjects) per cluster in group two. This value must be a positive number.

If you simply want a multiple of the value for group one, you would enter the multiple followed by M1, with no blanks. If you want to use M1 directly, you do not have to pre-multiply by 1. For example, all of the following are valid entries: 10 M1 2M1 0.5M1.

You can use a list of values such as 10 20 30 or M1 2M1 3M1.

Effect Size

Input Type

Indicate what type of values to enter to specify the differences. Regardless of the entry type chosen, the test statistics used in the power and sample size calculations are the same. This option is simply given for convenience in specifying the differences.

Effect Size - Differences

Only displayed if Input Type = "Differences"

D0 (Difference|H0 = P1.0 - P2)

Specify the difference between P1.0 and P2. This value is used with P2 to calculate the value of P1.0 using the formula: P1.0 = D0 + P2.

The power calculations assume that P1.0 is the value of the P1 under the null hypothesis. In this non-null case, the value of P1.0 is not equal to P2 as it is in the null case.

You may enter a range of values such as -.03 -.05 -.10 or -.05 to -.01 by .01. Differences must be between -1 and 1. D0 cannot take on the values -1, 0, or 1.

D1 (Difference|H1 = P1.1 - P2)

Specify the difference between P1.1 (P1 given the alternative hypothesis) and P2. This difference is used with P2 to calculate the value of P1 using the formula: P1.1 = D1 + P2.

The power calculations assume that P1.1 is the actual value of the proportion in group 1 (experimental or treatment group).

You may enter a range of values such as -.05 0 .5 or -.05 to .05 by .02. Actual differences must be between -1 and 1. They cannot take on the values -1 or 1.

Effect Size – Group 1 (Treatment)

Only displayed if Input Type = "Proportions"

P1.0 (Group 1 Proportion|H0)

Specify the value P1.0 which is the group 1 proportion given the null hypothesis.

The power calculations assume that P1.0 is the value of the P1 under the null hypothesis. In this non-null case, the value of P1.0 is not equal to P2 as it is in the null case.

Proportions must be between 0 and 1. They cannot take on the values 0 or 1. This value should not be set to exactly the value of P2.

You may enter a range of values such as 0.03 0.05 0.10 or 0.01 to 0.05 by 0.01.

P1.1 (Group 1 Proportion|H1)

Specify the value of P1.1 which is the group 1 proportion given the alternative hypothesis.

The power calculations assume that P1.1 is the actual value of the proportion in group 1 (experimental or treatment group).

Proportions must be between 0 and 1. They cannot take on the values 0 or 1.

You may enter a range of values such as 0.03 0.05 0.10 or 0.01 to 0.05 by 0.01.

Effect Size - Group 2 (Reference)

P2 (Group 2 Proportion)

Specify the value of p_2 , the control, baseline, or standard group's proportion. The null hypothesis is that the two proportions differ by a specified amount (See *Group 1 Proportion* below).

Since p_2 is a proportion, these values must be between 0 and 1.

You may enter a range of values such as 0.1 0.2 0.3 or 0.1 to 0.9 by 0.1.

Effect Size - Intracluster Correlation

ICC (Intracluster Correlation)

Enter a value (or range of values) for the intracluster correlation. This correlation may be thought of as the simple correlation between any two observations in the same cluster. It may also be thought of as the proportion of total variance in the observations that can be attributed to difference between clusters.

Although the actual range for this value is between 0 to 1, typical values range from 0.002 to 0.05.

Example 1 – Finding Power

A study is being designed to study the effectiveness of a new treatment. Historically, the standard treatment has enjoyed a 60% cure rate. The new treatment is more expensive and has more side effects than the standard treatment. Thus, the new treatment will be adopted if it is at least 8% better than the standard treatment.

The researchers will recruit patients from various hospitals. All patients at a particular hospital will receive the same treatment. They anticipate an average of 100 patients per hospital. Based on similar studies, they estimate the intracluster correlation to be 0.002.

The researchers plan to use the Farrington and Manning likelihood score test statistic to analyze the data. They want to study the power of the one-sided Farrington and Manning test at group cluster sizes ranging from 2 to 10 for detecting a difference of 0.08 when the actual cure rate of the new treatment ranges from 70% to 76%. The significance level will be 0.05.

Setup

This section presents the values of each of the parameters needed to run this example. First, from the PASS Home window, load the Non-Zero Null Tests for the Difference of Two Proportions in a Cluster-Randomized Design procedure window by expanding Proportions, then Two Proportions (Cluster-Randomized), then clicking on Test (Non-Zero Null), and then clicking on Non-Zero Null Tests for the Difference of Two Proportions in a Cluster-Randomized Design. You may then make the appropriate entries as listed below, or open Example 1 by going to the File menu and choosing Open Example Template.

<u>Option</u>	<u>Value</u>
Design Tab	
Solve For	Power
Alternative Hypothesis	One-Sided (H1: P1 - P2 > D0)
Test Type	Likelihood Score (Farr. & Mann.)
Alpha	0.05
K1 (Clusters in Group 1)	2 4 6 8 10
M1 (Items per Cluster in Group 1)	100
K2 (Clusters in Group 2)	K1
M2 (Items per Cluster in Group 2)	M1
Input Type	Differences
D0 (Difference H0 = P1.0 - P2)	0.08
D1 (Difference H1 = P1.1 - P2)	0.10 to 0.16 by 0.01
P2 (Group 2 Proportion)	0.6
ICC (Intracluster Correlation)	0.002

Output

Click the Calculate button to perform the calculations and generate the following output.

Numeric Results


Numeric Results Test Statistic: Likelihood Score Test (Farrington & Manning) Hypotheses: H0: P1 - P2 ≤ D0 vs. H1: P1 - P2 > D0 Group 1 Group 2 Intra-Clusters/ Clusters/ Group 2 Group 1 Group 1 Cluster Diff|H0 Diff|H1 **Items Items** Prop Prop|H0 Prop|H1 Corr. **Power** K1/M1 K2/M2 P2 P1.0 P1.1 D₀ **D1 ICC Alpha** 0.10341 2/100 2/100 0.6000 0.6800 0.7000 0.0800 0.1000 0.0020 0.050 0.13506 4/100 4/100 0.6000 0.6800 0.7000 0.0800 0.1000 0.0020 0.050 0.16346 6/100 6/100 0.6000 0.6800 0.7000 0.0800 0.1000 0.0020 0.050 0.19021 8/100 8/100 0.6000 0.6800 0.7000 0.0800 0.1000 0.0020 0.050 0.21588 10/100 10/100 0.6000 0.6800 0.7000 0.0800 0.1000 0.0020 0.050 0.14261 2/100 2/100 0.6000 0.6800 0.7100 0.0800 0.1100 0.0020 0.050 0.20377 4/100 4/100 0.6000 0.6800 0.7100 0.0800 0.1100 0.0020 0.050 0.25988 6/100 6/100 0.6000 0.6800 0.7100 0.0800 0.1100 0.0020 0.050 0.31271 8/100 8/100 0.6000 0.6800 0.7100 0.0800 0.1100 0.0020 0.050 0.36274 10/100 10/100 0.6000 0.6800 0.7100 0.0800 0.1100 0.0020 0.050 2/100 0.050 0.6800 0.0800 0.19129 2/100 0.6000 0.7200 0.1200 0.0020 0.29081 4/100 4/100 0.6000 0.6800 0.7200 0.0800 0.1200 0.0020 0.050 0.38063 6/100 0.6800 0.0800 0.050 6/100 0.6000 0.7200 0.1200 0.0020 0.46199 8/100 8/100 0.6000 0.6800 0.7200 0.0800 0.1200 0.0020 0.050 10/100 0.050 10/100 0.6800 0.0800 0.53513 0.6000 0.7200 0.1200 0.0020 0.24963 0.6000 0.6800 0.7300 0.0800 0.1300 2/100 2/100 0.0020 0.050 0.39323 4/100 4/100 0.6000 0.6800 0.7300 0.0800 0.1300 0.0020 0.050 0.1300 0.51582 6/100 6/100 0.6000 0.6800 0.0800 0.0020 0.050 0.7300 8/100 8/100 0.61848 0.6000 0.6800 0.7300 0.0800 0.1300 0.0020 0.050 0.70270 10/100 10/100 0.6000 0.6800 0.7300 0.0800 0.1300 0.050 0.0020 0.31705 2/100 2/100 0.6000 0.6800 0.7400 0.0800 0.1400 0.0020 0.050 0.50493 4/100 4/100 0.6000 0.6800 0.7400 0.0800 0.1400 0.0020 0.050 0.65065 6/100 6/100 0.6000 0.6800 0.7400 0.0800 0.1400 0.0020 0.050 0.75896 8/100 8/100 0.6000 0.6800 0.7400 0.0800 0.1400 0.0020 0.050 0.6800 0.7400 0.0800 0.83677 10/100 10/100 0.6000 0.1400 0.0020 0.050 **Report Definitions** H0 is an abbreviation for the null hypothesis. This is the hypothesis being evaluated by the statistical test. H1 is an abbreviation for the alternative hypothesis. This hypothesis gives the "true" parameter values. Power is the probability of rejecting a false null hypothesis. It should be close to one. K1 and K2 are the number of clusters in groups 1 and 2, respectively. M1 and M2 are the average number of items (subjects) per cluster in groups 1 and 2, respectively. P2 is the proportion for group 2, the standard, reference, baseline, or control group. P1.0 is the proportion for group 1 (the treatment group) assuming the null hypothesis (H0) P1.1 is the proportion for group 1 (the treatment group) assuming the alternative hypothesis (H1). D0 = P1.0 - P2 is the difference assuming the null hypothesis (H0). D1 = P1.1 - P2 is the actual difference at which the power is calculated. ICC is the intracluster correlation. Alpha is the probability of rejecting a true null hypothesis.

Summary Statements

Sample sizes of 200 in group 1 and 200 in group 2, which were obtained by sampling 2 clusters with 100 subjects each in group 1 and 2 clusters with 100 subjects each in group 2, achieve 10.341% power to detect a difference of 0.1000 when the null-hypothesized difference is 0.0800. The proportion in group 1 (the treatment group) is assumed to be 0.6800 under the null hypothesis and 0.7000 under the alternative hypothesis. The proportion in group 2 (the control group) is 0.6000. The test statistic used is the one-sided Likelihood Score Test (Farrington & Manning). The intracluster correlation is 0.0020, and the significance level of the test is 0.050.

This report shows the values of each of the parameters, one scenario per row. The total number of items sampled in group 1 is $N1 = K1 \times M1$. The total number of items sampled in group 2 is $N2 = K2 \times M2$.

Plots Section



The values from the table are displayed on the above charts. These charts give us a quick look at the sample sizes that will be required for various values of D1.

Example 2 – Finding the Sample Size (Number of Clusters)

Continuing with the scenario given in Example 1, the researchers want to determine the number of clusters necessary for each value of D1 when the target power is set to 0.80.

Setup

This section presents the values of each of the parameters needed to run this example. First, from the PASS Home window, load the Non-Zero Null Tests for the Difference of Two Proportions in a Cluster-Randomized Design procedure window by expanding Proportions, then Two Proportions (Cluster-Randomized), then clicking on Test (Non-Zero Null), and then clicking on Non-Zero Null Tests for the Difference of Two Proportions in a Cluster-Randomized Design. You may then make the appropriate entries as listed below, or open Example 2 by going to the File menu and choosing Open Example Template.

<u>Option</u>	<u>Value</u>
Design Tab	
Solve For	Sample Size (K1)
Alternative Hypothesis	One-Sided (H1: P1 - P2 > D0)
Test Type	Likelihood Score (Farr. & Mann.)
Power	0.80
Alpha	0.05
M1 (Items per Cluster in Group 1)	100
K2 (Clusters in Group 2)	K1
M2 (Items per Cluster in Group 2)	M1
Input Type	Differences
D0 (Difference H0 = P1.0 - P2)	0.08
D1 (Difference H1 = P1.1 - P2)	0.10 to 0.16 by 0.01
P2 (Group 2 Proportion)	0.6
ICC (Intracluster Correlation)	0.002

Output

Click the Calculate button to perform the calculations and generate the following output.

Numeric Results

Test Stati					& Manning	j)			
Hypotheses: H0: P1 - P2 ≤ D0 vs. H1: P1 - F Group 1 Group 2 Clusters/ Clusters/ Group 2 Group 1					Intra- Cluster				
Power	Items K1/M1	Items K2/M2	Prop P2	Prop H0 P1.0	Prop H1 P1.1	Diff H0 D0	Diff H1 D1	Corr. ICC	Alpha
0.80205	84/100	84/100	0.6000	0.6800	0.7000	0.0800	0.1000	0.0020	Alpha 0.050
0.80167	37/100	37/100	0.6000	0.6800	0.7100	0.0800	0.1100	0.0020	0.050
0.80764	21/100	21/100	0.6000	0.6800	0.7200	0.0800	0.1200	0.0020	0.050
0.82446	14/100	14/100	0.6000	0.6800	0.7300	0.0800	0.1300	0.0020	0.050
0.80122	9/100	9/100	0.6000	0.6800	0.7400	0.0800	0.1400	0.0020	0.050
0.82400	7/100	7/100	0.6000	0.6800	0.7500	0.0800	0.1500	0.0020	0.050
0.80374	5/100	5/100	0.6000	0.6800	0.7600	0.0800	0.1600	0.0020	0.050

The required sample size depends a great deal on the value of D1. The researchers should spend time determining the most appropriate value for D1.

Example 3 – Validation

This procedure uses the same mechanics as the Tests for Two Proportions in a Cluster-Randomized Design procedure. We refer the user to Example 4 of Chapter 230 for the validation.